

# C21orf62-AS1 Antibody



PACO35354

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

C21orf62-AS1

**Uniprot**

Q17RA5

**Synonyms:**

Putative uncharacterized protein C21orf62-AS1 (C21orf62 antisense RNA 1) (C21orf62 antisense gene protein 1), C21orf62-AS1, C21orf49

**Immunogen:**

Recombinant Human Putative uncharacterized protein C21orf62-AS1 protein (1-79AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A