

TCP10 Antibody



PACO35362

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TCP10

Uniprot

Q12799

Synonyms:

T-complex protein 10A homolog, TCP10, TCP10A

Immunogen:

Recombinant Human T-complex protein 10A homolog protein (134-353AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A