

ARL6IP1 Antibody, HRP conjugated



PACO35407

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in protein transport, membrane trafficking, or cell signaling during hematopoietic maturation.

Gene ID:

ARL6IP1

Uniprot

Q15041

Synonyms:

ADP-ribosylation factor-like protein 6-interacting protein 1 (ARL-6-interacting protein 1) (Aip-1) (Apoptotic regulator in the membrane of the endoplasmic reticulum), ARL6IP1, ARL6IP ARMER KIAA0069

Immunogen:

Recombinant Human ADP-ribosylation factor-like protein 6-interacting protein 1 protein (93-133AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A