

# RBMY1F Antibody, HRP conjugated



PACO35415

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

RNA-binding protein which may be involved in spermatogenesis. Required for sperm development, possibly by participating in pre-mRNA splicing in the testis.

**Gene ID:**

RBMY1F

**Uniprot**

Q15415

**Synonyms:**

RNA-binding motif protein, Y chromosome, family 1 member F/J (Y chromosome RNA recognition motif 2), RBMY1F; RBMY1J, YRRM2;

**Immunogen:**

Recombinant Human RNA-binding motif protein, Y chromosome, family 1 member F/J protein (1-496AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A