

LONRF1 Antibody



PACO35426

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

LONRF1

Uniprot

Q17RB8

Synonyms:

LON peptidase N-terminal domain and RING finger protein 1 (RING finger protein 191), LONRF1, RNF191


Immunogen:

Recombinant Human LON peptidase N-terminal domain and RING finger protein 1 protein (538-773AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human smooth muscle tissue using PACO35426 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO35426 at dilution of 1:100.