

# LONRF1 Antibody, HRP conjugated



PACO35427

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	LONRF1
<b>Source:</b>	<b>Uniprot</b>
Rabbit	Q17RB8
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	LON peptidase N-terminal domain and RING finger protein 1 (RING finger protein 191), LONRF1, RNF191
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human LON peptidase N-terminal domain and RING finger protein 1 protein (538-773AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A