

SMT3 Antibody



PACO35430

Product Information

Size:

50ug

Reactivity:

Saccharomyces cerevisiae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Not known; suppressor of MIF2 mutations.

Gene ID:

SMT3

Uniprot

Q12306

Synonyms:

Ubiquitin-like protein SMT3, SMT3

Immunogen:

Recombinant Saccharomyces cerevisiae Ubiquitin-like protein SMT3 protein (2-98AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A