

PACO35502

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

GTPase that associates with pre-60S ribosomal subunits in the nucleolus and is required for their nuclear export and maturation.

Gene ID:

GNL2

Uniprot

Q13823

Synonyms:

Nucleolar GTP-binding protein 2 (Autoantigen NGP-1), GNL2, NGP1

Immunogen:

Recombinant Human Nucleolar GTP-binding protein 2 protein (1-731AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A