

PACO35510

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Possible architectural role during spermatogenesis. May be involved in spermatid differentiation.

Gene ID:

CYLC2

Uniprot

Q14093

Synonyms:

Cylicin-2 (Cylicin II) (Multiple-band polypeptide II), CYLC2, CYL2

Immunogen:

Recombinant Human Cylicin-2 protein (1-156AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: CYLC2 antibody at 12 μ g/ml + 293T whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 40 kDa. Observed band size: 40 kDa.