

RHOH Antibody, FITC conjugated



PACO35528

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RHOH

Uniprot

Q15669

Synonyms:

Rho-related GTP-binding protein RhoH (GTP-binding protein TTF) (Translocation three four protein), RHOH, ARHH TTF

Immunogen:

Recombinant Human Rho-related GTP-binding protein RhoH protein (1-188AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A