

DGUOK Antibody, HRP conjugated



PACO35663

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Required for the phosphorylation of several deoxyribonucleosides and certain nucleoside analogs widely employed as antiviral and chemotherapeutic agents.

Gene ID:

DGUOK

Uniprot

Q16854

Synonyms:

Deoxyguanosine kinase, mitochondrial (dGK) (EC 2.7.1.113), DGUOK, DGK

Immunogen:

Recombinant Human Deoxyguanosine kinase, mitochondrial protein (40-277AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A