

# INPP5A Antibody, HRP conjugated



PACO35695

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Major isoenzyme hydrolyzing the calcium-mobilizing second messenger Ins(1,4,5)P<sub>3</sub>, this is a signal-terminating reaction.

**Gene ID:**

INPP5A

**Uniprot**

Q14642

**Synonyms:**

Type I inositol 1,4,5-trisphosphate 5-phosphatase (5PTase) (EC 3.1.3.56), INPP5A

**Immunogen:**

Recombinant Human Type I inositol 1,4,5-trisphosphate 5-phosphatase protein (1-409AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A