

TSC22D1 Antibody



PACO35718

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcriptional repressor. Acts on the C-type natriuretic peptide (CNP) promoter.

Gene ID:

TSC22D1

Uniprot

Q15714

Synonyms:

TSC22 domain family protein 1 (Cerebral protein 2) (Regulatory protein TSC-22) (TGFB-stimulated clone 22 homolog) (Transforming growth factor beta-1-induced transcript 4 protein), TSC22D1, KIAA1994 TGFB1I4 TSC22

Immunogen:

Recombinant Human TSC22 domain family protein 1 protein (1-144AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A