

# TSC22D1 Antibody, Biotin conjugated



PACO35721

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transcriptional repressor. Acts on the C-type natriuretic peptide (CNP) promoter.

**Gene ID:**

TSC22D1

**Uniprot**

Q15714

**Synonyms:**

TSC22 domain family protein 1 (Cerebral protein 2) (Regulatory protein TSC-22) (TGFB-stimulated clone 22 homolog) (Transforming growth factor beta-1-induced transcript 4 protein), TSC22D1, KIAA1994 TGFB1I4 TSC22

**Immunogen:**

Recombinant Human TSC22 domain family protein 1 protein (1-144AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A