

FAM187B Antibody



PACO35734

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

FAM187B

Uniprot

Q17R55

Synonyms:

Protein FAM187B (Transmembrane protein 162), FAM187B, TMEM162


Immunogen:

Recombinant Human Protein FAM187B protein (18-335AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO35734 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO35734 at dilution of 1:100.