

MFAP5 Antibody, FITC conjugated



PACO35752

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of the elastin-associated microfibrils.

Gene ID:

MFAP5

Uniprot

Q13361

Synonyms:

Microfibrillar-associated protein 5 (MFAP-5) (MP25) (Microfibril-associated glycoprotein 2) (MAGP-2), MFAP5, MAGP2

Immunogen:

Recombinant Human Microfibrillar-associated protein 5 protein (22-173AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A