

DIP2A Antibody



PACO35770

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May provide positional cues for axon pathfinding and patterning in the central nervous system.

Gene ID:

DIP2A

Uniprot

Q14689

Synonyms:

Disco-interacting protein 2 homolog A (DIP2 homolog A), DIP2A, C21orf106 DIP2 KIAA0184


Immunogen:

Recombinant Human Disco-interacting protein 2 homolog A protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO35770 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO35770 at dilution of 1:100.