

ARHGAP19 Antibody



PACO35782

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.

Gene ID:

ARHGAP19

Uniprot

Q14CB8

Synonyms:

Rho GTPase-activating protein 19 (Rho-type GTPase-activating protein 19), ARHGAP19


Immunogen:

Recombinant Human Rho GTPase-activating protein 19 protein (1-494AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO35782 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO35782 at dilution of 1:100.