

NGF Antibody, HRP conjugated



PACO35851

Product Information

Size:

50ug

Reactivity:

Pig

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Nerve growth factor is important for the development and maintenance of the sympathetic and sensory nervous systems. Extracellular ligand for the NTRK1 and NGFR receptors, activates cellular signaling cascades through those receptor tyrosine kinase to regulate neuronal proliferation, differentiation and survival. Inhibits metalloproteinase dependent proteolysis of platelet glycoprotein VI.

Gene ID:

NGF

Uniprot

Q29074

Synonyms:

Beta-nerve growth factor (Beta-NGF) (Fragment), NGF, NGFB

Immunogen:

Recombinant Pig Beta-nerve growth factor protein (110-229AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A