

PRICKLE4 Antibody



PACO35886

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

PRICKLE4

Uniprot

Q2TBC4

Synonyms:

Prickle-like protein 4 (Overexpressed breast tumor protein), PRICKLE4, C6orf49 OEBT


Immunogen:

Recombinant Human Prickle-like protein 4 protein (1-240AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO35886 at dilution of 1:100.