

PRICKLE4 Antibody, Biotin conjugated



PACO35889

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PRICKLE4

Uniprot

Q2TBC4

Synonyms:

Prickle-like protein 4 (Overexpressed breast tumor protein), PRICKLE4, C6orf49 OEBT

Immunogen:

Recombinant Human Prickle-like protein 4 protein (1-240AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A