

C2orf76 Antibody



PACO35914

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

C2orf76

Uniprot

Q3KRA6

Synonyms:

UPF0538 protein C2orf76, C2orf76


Immunogen:

Recombinant Human UPF0538 protein C2orf76 protein (1-126AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO35914 at dilution of 1:100.