

ALDH1L2 Antibody, HRP conjugated



PACO35923

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ALDH1L2

Uniprot

Q3SY69

Synonyms:

Mitochondrial 10-formyltetrahydrofolate dehydrogenase (Mitochondrial 10-FTHFDH) (mtFDH) (EC 1.5.1.6) (Aldehyde dehydrogenase family 1 member L2), ALDH1L2

Immunogen:

Recombinant Human Mitochondrial 10-formyltetrahydrofolate dehydrogenase protein (263-562AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A