

TSGA10IP Antibody



PACO35930

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TSGA10IP

Uniprot

Q3SY00

Synonyms:

Testis-specific protein 10-interacting protein (Tsga10-interacting protein), TSGA10IP

Immunogen:

Recombinant Human Testis-specific protein 10-interacting protein (6-302AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A