

# SLC10A6 Antibody



PACO35970

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transports sulfoconjugated steroid hormones, as well as tauro lithocholic acid, 3-sulfate and sulfoconjugated pyrenes in a sodium-dependent manner.

**Gene ID:**

SLC10A6

**Uniprot**

Q3KNW5

**Synonyms:**

Solute carrier family 10 member 6 (Sodium-dependent organic anion transporter), SLC10A6, SOAT

**Immunogen:**

Recombinant Human Solute carrier family 10 member 6 protein (312-377AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A