

HYDIN Antibody



PACO36006

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Required for ciliary motility.

Gene ID:

HYDIN

Uniprot

Q4G0P3

Synonyms:

Hydrocephalus-inducing protein homolog, HYDIN, HYDIN1 KIAA1864

Immunogen:

Recombinant Human Hydrocephalus-inducing protein homolog protein (101-400AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A