

# FAM153C Antibody



PACO36038


---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	FAM153C
<b>Source:</b>	<b>Uniprot</b>
Rabbit	Q494X1
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Protein FAM153C, FAM153C
<b>Applications:</b>	<b>Immunogen:</b>
ELISA, IHC	Recombinant Human Protein FAM153C protein (1-144AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
ELISA:1:2000-1:10000, IHC:1:20-1:200	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO36038 at dilution of 1:100.