

GTDC1 Antibody



PACO36050

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

GTDC1

Uniprot

Q4AE62

Synonyms:

Glycosyltransferase-like domain-containing protein 1 (Mannosyltransferase-like protein Xa) (Mat-Xa), GTDC1


Immunogen:

Recombinant Human Glycosyltransferase-like domain-containing protein 1 protein (1-292AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: GTDC1 antibody at 4 μ g/ml + A549 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 53, 44, 34, 50, 38 kDa. Observed band size: 53 kDa.