

FAM135B Antibody



PACO36058

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

FAM135B

Uniprot

Q49AJ0

Synonyms:

Protein FAM135B, FAM135B, C8orfK32


Immunogen:

Recombinant Human Protein FAM135B protein (801-1100AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO36058 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO36058 at dilution of 1:100.