

VWA3B Antibody



PACO36078

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

VWA3B

Uniprot

Q502W6

Synonyms:

von Willebrand factor A domain-containing protein 3B (VWA domain-containing protein 3B), VWA3B


Immunogen:

Recombinant Human von Willebrand factor A domain-containing protein 3B protein (1-160AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human skin tissue using PACO36078 at dilution of 1:100.