

SLC44A4 Antibody, Biotin conjugated



PACO36093

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SLC44A4

Uniprot

Q53GD3

Synonyms:

Choline transporter-like protein 4 (Solute carrier family 44 member 4) (Thiamine pyrophosphate transporter 1) (hTPPT1), SLC44A4, C6orf29 CTL4 NG22 TPPT1

Immunogen:

Recombinant Human Choline transporter-like protein 4 protein (56-229AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A