

C2CD6 Antibody



PACO36098

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

C2CD6

Uniprot

Q53TS8

Synonyms:

C2 calcium-dependent domain-containing protein 6 (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 11 protein), C2CD6, ALS2CR11


Immunogen:

Recombinant Human C2 calcium-dependent domain-containing protein 6 protein (401-623AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO36098 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human melanoma using PACO36098 at dilution of 1:100.