

ssaA1 Antibody



PACO36102

Product Information

Size:

50ug

Reactivity:

Staphylococcus aureus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Not known; immunogenic protein.

Gene ID:

ssaA1

Uniprot

Q5HCY4

Synonyms:

Staphylococcal secretory antigen ssaA1, ssaA1

Immunogen:

Recombinant Staphylococcus aureus Staphylococcal secretory antigen ssaA1 protein (27-255AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A