

fabD Antibody, HRP conjugated



PACO36131

Product Information

Size:

50ug

Reactivity:

Staphylococcus aureus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

fabD

Uniprot

Q5HGK3

Synonyms:

Malonyl CoA-acyl carrier protein transacylase (MCT) (EC 2.3.1.39), fabD

Immunogen:

Recombinant Staphylococcus aureus Malonyl CoA-acyl carrier protein transacylase protein (1-308AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A