

TTC32 Antibody




PACO36134


Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	TTC32
Source:	Uniprot
Rabbit	Q5I0X7
Isotype:	Synonyms:
IgG	Tetratricopeptide repeat protein 32 (TPR repeat protein 32), TTC32
Applications:	Immunogen:
ELISA, IHC	Recombinant Human Tetratricopeptide repeat protein 32 protein (1-151AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, IHC:1:20-1:200	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO36134 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO36134 at dilution of 1:100.