

# CFAP157 Antibody



PACO36138

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

CFAP157

**Uniprot**

Q5JU67

**Synonyms:**

Cilia- and flagella-associated protein 157, CFAP157, C9orf117

**Immunogen:**


Recombinant Human Cilia- and flagella-associated protein 157 protein (1-200AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human lung tissue using PACO36138 at dilution of 1:100.