

C4orf46 Antibody



PACO36150

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

C4orf46

Uniprot

Q504U0

Synonyms:

Renal cancer differentiation gene 1 protein, C4orf46, RCDG1


Immunogen:

Recombinant Human Renal cancer differentiation gene 1 protein (1-113AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO36150 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO36150 at dilution of 1:100.