

# RFTN2 Antibody, HRP conjugated



PACO36155

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Upon bacterial lipopolysaccharide stimulation, mediates clathrin-dependent internalization of TLR4 in dendritic cells, resulting in activation of TICAM1-mediated signaling and subsequent IFNB1 production. May regulate B-cell antigen receptor-mediated signaling.

**Gene ID:**

RFTN2

**Uniprot**

Q52LD8

**Synonyms:**

Raftlin-2 (Raft-linking protein 2), RFTN2, C2orf11

**Immunogen:**

Recombinant Human Raftlin-2 protein (1-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A