

UNC50 Antibody, HRP conjugated



PACO36163

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in cell surface expression of neuronal nicotinic receptors. Binds RNA.

Gene ID:

UNC50

Uniprot

Q53HI1

Synonyms:

Protein unc-50 homolog (Periodontal ligament-specific protein 22) (PDLs22) (Protein GMH1 homolog) (hGMH1) (Uncoordinated-like protein), UNC50, UNCL

Immunogen:

Recombinant Human Protein unc-50 homolog protein (15-82AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A