

PACO36182

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

GCSAML

Uniprot

Q5JQS6

Synonyms:

Germinal center-associated signaling and motility-like protein, GCSAML, C1orf150


Immunogen:

Recombinant Human Germinal center-associated signaling and motility-like protein (1-135AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO36182 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO36182 at dilution of 1:100.