

PACO36202

Product Information

Size:

50ug

Reactivity:

Staphylococcus aureus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the condensation reaction of fatty acid, synthesis by the addition to an acyl acceptor of two carbons from malonyl-ACP.

Gene ID:

fabF

Uniprot

Q5HHA1

Synonyms:

3-oxoacyl-[acyl-carrier-protein] synthase 2 (EC 2.3.1.179) (3-oxoacyl-[acyl-carrier-protein] synthase II) (Beta-ketoacyl-ACP synthase II) (KAS II), fabF

Immunogen:

Recombinant Staphylococcus aureus 3-oxoacyl-[acyl-carrier-protein] synthase 2 protein (1-414AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A