

# GLE1 Antibody, FITC conjugated



PACO36212

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Required for the export of mRNAs containing poly(A) tails from the nucleus into the cytoplasm. May be involved in the terminal step of the mRNA transport through the nuclear pore complex (NPC).

**Gene ID:**

GLE1

**Uniprot**

Q53GS7

**Synonyms:**

Nucleoporin GLE1 (hGLE1) (GLE1-like protein), GLE1, GLE1L

**Immunogen:**

Recombinant Human Nucleoporin GLE1 protein (141-380AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A