

choD Antibody



PACO36226

Product Information

Size:

50ug

Reactivity:

Mycobacterium tuberculosis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the oxidation and isomerization of cholesterol to cholestenone (4-cholesten-3-one), which is an initial step in the cholesterol degradation process. Required for virulence.

Gene ID:

choD

Uniprot

P9WMV8

Synonyms:

Cholesterol oxidase (EC 1.1.3.6) (Cholesterol isomerase) (EC 5.3.3.1), choD

Immunogen:

Recombinant Mycobacterium tuberculosis Cholesterol oxidase protein (1-578AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A