

BEX5 Antibody



PACO36234

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

extracellular exosome.

Gene ID:

BEX5

Uniprot

Q5H9J7

Synonyms:

Protein BEX5 (Brain-expressed X-linked protein 5) (NGFRAP1-like protein 1) (Nerve growth factor receptor-associated protein 2), BEX5, NADE2 NGFRAP1L1


Immunogen:

Recombinant Human Protein BEX5 protein (1-111AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO36234 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO36234 at dilution of 1:100.