

MAP3K19 Antibody



PACO36254

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MAP3K19

Uniprot

Q56UN5

Synonyms:

Mitogen-activated protein kinase kinase kinase 19 (EC 2.7.11.1) (Regulated in COPD, protein kinase) (SPS1/STE20-related protein kinase YSK4), MAP3K19, RCK YSK4

Immunogen:

Recombinant Human Mitogen-activated protein kinase kinase kinase 19 protein (45-168AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A