

DALRD3 Antibody



PACO36258

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

DALRD3

Uniprot

Q5D0E6

Synonyms:

DALR anticodon-binding domain-containing protein 3, DALRD3

Immunogen:

Recombinant Human DALR anticodon-binding domain-containing protein 3 protein (182-481AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A