

PAGE2B Antibody



PACO36270

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

PAGE2B

Uniprot

Q5JRK9

Synonyms:

Putative G antigen family E member 3 (Prostate-associated gene 2B protein) (PAGE-2B), PAGE2B, GAGEE3


Immunogen:

Recombinant Human Putative G antigen family E member 3 protein (1-111AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO36270 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO36270 at dilution of 1:100.