

TMTC4 Antibody



PACO36326

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

TMTC4

Uniprot

Q5T4D3

Synonyms:

Transmembrane and TPR repeat-containing protein 4, TMTC4


Immunogen:

Recombinant Human Transmembrane and TPR repeat-containing protein 4 protein (245-376AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO36326 at dilution of 1:100.