

ATAD3C Antibody



PACO36362

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

ATAD3C

Uniprot

Q5T2N8

Synonyms:

ATPase family AAA domain-containing protein 3C, ATAD3C


Immunogen:

Recombinant Human ATPase family AAA domain-containing protein 3C protein (1-122AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO36362 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO36362 at dilution of 1:100.