

ASPHD1 Antibody



PACO36370

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

ASPHD1

Uniprot

Q5U4P2

Synonyms:

Aspartate beta-hydroxylase domain-containing protein 1 (EC 1.14.11. -), ASPHD1


Immunogen:

Recombinant Human Aspartate beta-hydroxylase domain-containing protein 1 protein (96-390AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO36370 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skin tissue using PACO36370 at dilution of 1:100.